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CertoClav CV-EL

12/18 L

Laborautoklav

(User manual English: page 15)

Inhalt

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(German version: page 1)

Content

Content	16
Introduction.....	17
Declaration of conformity CV-EL.....	18
General safety instructions.....	19
Autoclave – Safegurads	20
Technical Data	21
<i>Overview – Autoclave</i>	<i>21</i>
<i>Directives.....</i>	<i>21</i>
<i>Lables</i>	<i>22</i>
Scope of Delivery.....	23
Setup	23
<i>Autoclave.....</i>	<i>23</i>
Operation	24
<i>Before first use.....</i>	<i>24</i>
<i>Loading the autoclave.....</i>	<i>24</i>
<i>Close and Locking.....</i>	<i>24</i>
<i>Sterilization / Autoclaving.....</i>	<i>25</i>
<i>End of cycle – Cooling down.....</i>	<i>25</i>
<i>Opening</i>	<i>25</i>
TOTAL BATCH TIME	26
Autoclave - cleaning and maintenance.....	26
<i>Replacement of the lid gasket.....</i>	<i>27</i>
Spare parts	27
Electric schematic.....	28

Introduction

Congratulations on the purchase of the CertoClav CV-EL.

The CV-EL is a classic among our laboratory autoclaves. It impresses with its sturdy construction and simple operation. It belongs to the class N steam sterilizers. The venting is done by a flow process.

The device is controlled by a valve. This valve is available in two versions, each with two temperature levels. You can choose between 115/121°C or 125/140°C. The temperature is set with the control valve in the middle of the lid. The valves cannot be interchanged, as the thermostat of the unit is harmonized with and calibrated to the respective valve!

The device has been manufactured using modern manufacturing methods and has been 100% tested by our quality control department. The device and all parts, not subject to wear and tear, are covered by a warranty, for 2 years, from the date of purchase. Deposits and discolorations result from the test runs carried out. No warranty is given in case of improper use, e.g. if the charge contains corrosive material.

The CertoClav CV-EL is a technical device and requires certain expertise in operation to achieve error-free results. Please read this manual to learn about all possibilities and handling.

Please follow the safety warnings.

Should you have any questions, please do not hesitate to contact us.

By e-mail: office@certoclav.com or phone: +43 732 674 278-0.

Your CertoClav Team

Declaration of conformity CV-EL

This is to certify, that the **CERTOCLAV** Product Range CV-EL complies with the following regulations of the European Community:

2014/35/EU Low Voltage Directive

Associated Standards:

EN 61010-2-040:2015 Safety requirements for electrical equipment for measurements, control and laboratory Part 2-040: Particular requirements for sterilizers and washer-disinfectors used to treat medical materials.

EN 61010-1:2010 Safety requirements for electrical equipment for measurements, control and laboratory use – Part 1: General requirements

The CertoClav pressure equipment is constructed according Low Voltage Directive. 2014/30/EU Electromagnetic compatibility

Associated Standards:

EN 61326-1:2013 Electrical Equipment for measurements, control and laboratory use - EMC Requirements- Part 1: General requirements (provided replacement: EDIN EN 61326-1:2018)

EN 61000-3-2:2014 Electromagnetic Compatibility (EMC) — Part 3- 2: Limits —

Limits for harmonic current emissions (equipment input current equal to or less than 16 A per phase)

EN 61000-3-3:2013 Electromagnetic compatibility — Part 3- 3: Limits — Limits for voltage fluctuation and flicker on the public low-voltage system (equipment having an input current equal to or less 16 A per phase) and not subject to conditional connection.

2011/65/EU Restriction of the use of certain hazardous substances (RoHS) Directive

Associated Standards:

EN IEC 63000:2018 Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances

Manufacturing quality:

CertoClav manufactures according to modern manufacturing methods using high-quality materials. Decades of experience and experienced employees guarantee high product quality and maximum customer service. All devices are 100% tested for function and quality. CertoClav products are manufactured exclusively in Austria.

Manufacturer:

CERTOCLAV Sterilizer GmbH
A-4060 Leonding / Austria

Michael Dirix CEO

Leonding, June 2022

General safety instructions

-  Only use the devices in accordance with the user manual. Keep the user manual in a place accessible to all users of the device.
-  Never tamper with safety system and safeguards. Do NOT use the lid and chamber of different series together.
-  The surface of the device (metal and plastic) becomes hot during operation. Use suitable protective gloves.
-  This device must NOT be installed in an area, that is subject to the danger of explosions.
-  Only use the device with the intended mains voltage.
-  Only intended for use in laboratories by technically trained personnel. Do NOT use in household settings.
-  Never operate the unit without water.
-  Fill level for 12 liter unit: min. 1 liter (max. 1.5l); for 18 liter unit: min. 1 liter (max. 2l).
-  Do NOT autoclave caustic or corrosive charges or saline solutions. This will cause corrosion of the unit, which is made of an aluminum alloy.
-  Demineralized water can cause pitting corrosion (mix with tap water to avoid this).
-  When not in use, be sure to remove the water from the chamber (to avoid deposits and pitting corrosion!).
-  Do NOT heat any tightly closed or sealed vessels (especially glass bottles), as there is no pressure compensation.
-  Use only water to generate steam. We recommend distilled water. DO NOT use demineralized water from an ion exchanger, this will cause pitting corrosion. If only demineralized water is available, mix with tap water (approx. 1/3)!
-  DO NOT release steam in the direction of people - risk of burning. Fit a temperature-resistant hose (ø inside 10 mm) to the exhaust steam spigot, so that the hot steam can be discharged directly and safely into a sink/ container. Article no. 8582016 - Hose 1.0m.
-  When autoclaving liquids, there is the possibility of boiling the liquid (boiling delay), if the steam is quickly released manually. To avoid this, the autoclave should be cooled down to 70°C at room temperature without opening a valve.

-  Do not transport the device when it is filled.
-  The device may only be repaired or maintained by CertoClav Sterilizer GmbH or by eligible personnel, authorized by CertoClav Sterilizer GmbH (specialized company or customer service), using original spare parts only and following the service instructions.
-  Avoid contact of the aluminum chamber with stainless steel due to chemical processes.

Autoclave – Safeguards

Secure locking mechanism and steam release valve within the lid

To prevent opening under pressure and pressure build-up, if the device is not properly closed. Turning the orange knob to the "OPEN" position automatically opens the steam release valve and releases steam. As long as the device is not properly locked, the steam release valve remains open and prevents a pressure build-up.

Pressure Control valve:

If the maximum working pressure of 2.7 bar is exceeded by max. 10%, this pressure control valve opens.

Safety valve:

If the pressure control valve fails, the safety valve opens at a pressure of 3.6 to 4 bar (and closes again after reaching the operating pressure). The safety valve must then be replaced in any case.

Overpressure safety recess

In the chamber rim, left to the device handle, there is a notch. If all other safeguards should fail and the pressure reaches around 4,7 bar, the seal inside the lid will be pressed through this overpressure safety notch with a loud noise. This way the device is no longer sealed, the pressure and temperature will start to sink.

If this occurs do NOT operate the unit anymore and contact your dealer or the manufacturer!

-  If this occurs do not operate the unit anymore and contact your dealer or the manufacturer for examination and check. It is dangerous to work with such a unit.

Technical Data

Overview – Autoclave

	CV-EL 12L	CV-EL 18L
Ambient temperature	-10°...+40°C	
Chamber volume	12 Liter	18 Liter
Temperature / Pressure	140°C; 2,7 bar / 125°C; 1,4 bar / 121°C; 1,2 bar / 115°C; 0,7 bar	
Chamber measurements	Ø 24cm, 26 cm high	Ø 24cm, 39 cm high
Diagonal	35 cm	46,5 cm
Max. Working pressure	2,7 bar	
Test pressure	4,1 bar	
Weight	7 kg	9 kg
Material	Aluminum No. 3103 according to DIN1725	
Voltage	230V ~10% 50 Hz	
Power	Heating: 1900 W	
Pressure equipment designed according to	Low-Voltage Directive 2014/35/EU	

Directives

- **EG-Low-voltage-directive 2014/35/EU**

CertoClav CV-EL is constructed according to the **Low-Voltage-Directive 2014/35/EU**.

Tested in accordance with:

EN61010-1:2001 (IEC 1010-1:2001) Safety requirements for electrical equipment for measurement, control and laboratory use – Part 1: General requirements

EN61010-2-040:2005 (IEC 1010-2-040:2005) Particular requirements for sterilizers and washer-disinfectors used to treat medical materials

- **Electromagnetic Compatibility Directive 2014/30/EU**

Associated standards:

EN 61326-1:2013 Electrical Equipment for measurements, control and laboratory use - EMC-Requirements Part 1: General requirements

EN 61000-3-2:2006 + A1:2009 + A2:2009 Electromagnetic Compatibility (EMC) — Part 3- 2: Limits - Limits for harmonic current emissions (Equipment input current equal to or less than ≤ 16A per phase)

EN 61000-3-3:2013 Electromagnetic Compatibility — Part 3- 3: Limits — Limits for voltage fluctuation and flicker on the public low-voltage system (equipment having an input current to or less ≤ 16 A per phase), and not subject to conditional connection.

EN 61010-1:2010


- **RoHS 2011/65/EU: Directive on the restriction of the use of certain hazardous substances**

Expert opinion on germicidal effectiveness and deaeration performance prepared by the Austrian Federal State Bacteriological-Serological Research Institute in Vienna, Institute for Applied Hygiene in Graz and ANAWA in Munich.

CERTOCLAV CV-EL is **not** suitable for sterilization of tightly sealed bottles and loads with aggressive media. For porous loads the **CERTOCLAV** CV-EL CANNOT be recommended. Smaller amounts of textiles with low density can be sterilized.

We recommend the use of sterilization indicators.

Lables

	EAG Directive 2012/19/EU: Old electrical equipment is taken back by the manufacturer or can be handed in at local collection points
	ATTENTION: Always follow the user manual
	ATTENTION: Hot surface! Do not touch!
	The CE mark certifies the conformity of the device to the directives NSR 2014/35 EU + EMC 2014/30 EU, and EN 61010-1:2010.
	Short time operation: 30 min

Scope of Delivery

After opening the packaging, please inspect the product for eventual transport damage. In case of damaged items, contact your dealer and do NOT operate the unit for safety reasons. Please also check whether the full scope of delivery is included.

The following parts are delivered with the autoclave, please check:

- Autoclave chamber CV-EL and lid
- Instrument tray with water level indicator
- Manometer with gasket and nut
- Steam release spigot with gasket and nut
- Power cable
- User manual and warranty certificate

If anything listed above is missing, please contact CertoClav Sterilizer GmbH.

Additional we offer the following accessories:

- Art. no. 8583010 Wire basket Ø 23, stainless steel
- Art. no. 8582020 Condensation collection container (5L) and hose (1,0m)
- Art. no. 8582016 Hose (1,0m) for mounting to steam release cock (Ø 10mm)
- Art. no. 1250499 Service Set 1 with all seals
- Art. no. 8500482 Water distiller. We recommend using distilled water

For further information about our product range visit: www.certoclav.com

Setup

Autoclave

- The installation is straight forward, either on a table or on the floor, so that an ideal working height is given.
- Unpack the manometer and the steam release spigot, insert them with the seal through the opening on the left and right side of the lid and fix them with the nut from the inside of the lid.
- We recommend mounting a temperature resistant hose to the steam release spigot (Art. No. 8582016), inner diameter 10 mm.
- Plug the power cable into the socket on the device base. The ON/OFF switch is located on the front.

Operation

Before first use

-  READ INSTRUCTIONS BEFORE FIRST USE!
-  The first cycle must be run without load and water only!

Loading the autoclave

- Fill the chamber with 1 to 1.5 liters of distilled water. (if demineralized, mix with tap water)
ATTENTION: No demineralized water from an ion exchanger, it causes pitting. If only this is available, mix with tap water (approx. 1/3)!
 - Place the instrument plate, container or basket in the chamber and insert the load. It is important that the steam can circulate freely.
 - ATTENTION: No stainless-steel parts may come into direct contact with the chamber - causes pitting!
-  **Do not autoclave tightly sealed containers (especially glass bottles) since there is no pressure compensation.**

Close and Locking

- Turn the orange knob on the lid handle to the "OPEN" position.
 - Open the steam release spigot (turn the black knob counterclockwise).
 - Connect the hose to prevent danger from hot steam or water.
 - Place the lid on the chamber in such a way that the arrows on the lid and in the middle of the handles are aligned. Then turn the lid clockwise until the lid and chamber handles are aligned.
 - Turn the orange knob to "CLOSE" (= lock).
 - Set the red drag pointer on the manometer to "0".
 - Adjust the pressure control valve to the desired temperature by turning the valve cap as far as it will go. Depending on the valve: for 125°C / 115°C turn counterclockwise, for 140° C / 121°C turn clockwise. Never open the small screw!
 - Adjust the temperature at the regulator knob according to the valve.
 - Switch on the device with the green on/off switch.
-  Long-term cycles of more than two hours may result in dry boiling and overheating. This will damage the autoclave.

Sterilizing / Autoclaving

Heating up

6 to 8 minutes after switching on, steam begins to escape from the open steam release spigot. To remove the air effectively, vent for 3 minutes, then close the steam release spigot (turn the knob clockwise).

CAUTION: Without hose connection danger from hot steam!

After approx. 3 minutes the pressure and desired temperature are reached.

The pressure is controlled via the central pressure control valve, steam is blown off regularly. In CV-EL units the energy supply is controlled via a thermostat. To enable interaction, set the same temperature for the pressure control valve (115/121°C or 125/140°C) and the thermostat.

Sterilization / Holding time

The start is indicated by the first escape of steam from the pressure control valve. The set temperature should be reached and the pressure indicator on the manometer should be in the red or blue marking. If necessary, set a timer.

After the sterilization time has expired, switch off with the green on/off switch.

ATTENTION: Danger of dry boiling and overheating if sterilization time is too long. The autoclave will be damaged.

End of cycle – Cooling down

Make sure that the pointer of the manometer points at zero, check the thermometer and the indicator label on the chamber.

When heating glass bottles, before opening the autoclave, wait until the temperature has fallen below 80°C, for media or liquids below 70°C (after about 40 minutes). Otherwise, there is a risk of boiling delay. After the temperature has fallen, open the steam release spigot.

Opening

Open the steam release spigot: Turn the orange knob to "OPEN", this will open the exhaust valve, steam will escape. Turn the lid counterclockwise and remove it. To improve drying, remove the basket from the chamber and allow it to cool down to room temperature before removing the load.

INCOMPLETE CYCLE: If, for any reason, the cycle is interrupted (by switching off the device or power failure), the sterilization process has not been carried out correctly and must be repeated. Before opening, follow the safety instructions for cooling down.

TOTAL BATCH TIME

A complete working cycle consists of:

Heating time:	depending on the load and size of the unit	6 - 8 minutes
Venting time:	minimum	3 minutes
Pressure build-up:	about	3 minutes

Sterilization/holding time and temperature: Depending on application and load

Follow the manufacturer's specifications e.g. 140°C - 5 min, 125°C - 15 min, 121°C - min. 15-20 minutes

Cooling time and drying: full bodies can be removed after pressure reduction

Glass at 80°C: after about 30 minutes

Liquids at 70°C: after about 40 minutes.

Autoclave - cleaning and maintenance

For cleaning, simply wipe inside and outside with a damp cloth. Do NOT immerse the lid and chamber in water. Before cleaning/disinfecting, pull the power plug and clean/disinfect the unit only when cold.

Disinfectants with the following ingredients are suitable for disinfecting the surface of the housing:

Ethanol, propanol, isopropanol, glutaraldehyde, quaternary ammonium compounds. It is essential to follow the application instructions of the manufacturers of cleaning and disinfecting agents.

Deposits and discolorations are caused by normal water. These deposits can clog the pressure control valve. Therefore, always use distilled water (if demineralized, mix with approx. 1/3 tap water).

To remove discoloration, add some citric acid to the water in the chamber and boil for about 10 minutes without putting on the lid. Then empty immediately.


Only use original spare parts for maintenance! Repairs may only be carried out by the manufacturer or by expressly authorized personnel.

It is recommended to have the appliance inspected at an interval of 2 years by your dealer or the manufacturer.

When the unit is not in use, the water must be removed from the chamber to avoid deposits and pitting corrosion!

Replacement of the lid gasket


1. Remove the old gasket and clean the groove.
2. Oil the groove and the new gasket with cooking oil
3. Press gasket into groove ribbed side pointing downward (letters show to the middle of the lid) Do not stretch the gasket when inserting, as it will get too long. Rather try to squeeze it a little bit, it will seal better.

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
Art. No. 1250499 Service Set 1 - The Service Set 1 contains all gaskets and valves which are necessary for the maintenance of your autoclave lid. Seals must be replaced regularly. **During long-term cycles the wear of seals is considerably higher.** Replacement may be necessary here after only 6 months.

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