

Alcian blue pH 1,0 Mowry technical's information

Technical's card code 09-106

Product code 09-106

Pack 100 ml or on request

Stability of product properly conserved in a dark place at 15-25°C 24 months.




#### Suppliers

Company identification.

DDKItalia S.r.l

Via Marche 19-27029 Vigevano (I)

[info@ddkitalia.com](mailto:info@ddkitalia.com) • [www.ddkitalia.com](http://www.ddkitalia.com)

in case of emergency UE number		112
in case of emergency UK number		999
en cas d'urgence Suisse number		145

#### Principle

Combination of PAS staining with Alcian blue allows the additional identification of acid mucosubstances (glycosaminoglycans). For PAS staining the material is treated with periodic acid during which process the 1,2-glycols are oxidized to aldehydes groups. With Schiff's reagent, aldehydes give a bright red color reaction. With non-substituted polysaccharides, neutral mucopolysaccharides, mucoproteins and glycoprotein, glycolipids and phospholipids the PAS staining gives a specific color reaction. Acid mucosubstances can be demonstrated by the non-specific method with nuclear fast red counterstaining, without differentiation of carboxyl and sulphate groups.

#### Sample material

Sections of formalin-fixed, paraffin-embedded tissue or cell smears are used as starting material.

3-5 µm thick paraffin sections

#### Method

Alcian blue staining for acid mucosubstances

Deparaffinise the sections in the typical manner and rehydrate

Alcian blue solution	5 minute, Running tap water	3 minute
----------------------	-----------------------------	----------

Rinse in distilled water

Nuclear fast red solution	10 minute. Running tap water	3 minute
---------------------------	------------------------------	----------

Rinse in distilled water

Ascending alcohol series

Clear in xylene

Mount with DdMount

#### Result

Nuclei	red
Acid mucosubstances	light blue

#### Method

Procedure Alcian blue-PAS staining

Deparaffinize the sections in the typical manner and rehydrate

Alcian blue solution	5 minute. Running tap water	3 minute
----------------------	-----------------------------	----------

Rinse in distilled water

Periodic acid	10 minute. Running tap water	3 minute
---------------	------------------------------	----------

Rinse in distilled water

Schiff's reagent	15 minute
------------------	-----------

Running tap water	3 minute.
-------------------	-----------

Rinse in distilled water

Gill III haematoxylin solution	20 second. Running tap water	3 minute
--------------------------------	------------------------------	----------

Ascending alcohol series

Clear in xylene

Mount with DdMount

Alcian blue pH 1,0 Mowry technical's information  
 Technical's card code 09-106  
 Product code 09-106

### Result

Nuclei	blue
Acid mucosubstances	light blue
Polysaccharides, neutral mucopolysaccharides	purple

### Sample preparation.

All samples should be treated using state-of-the-art technology. All samples must be clearly labelled. Suitable instruments must be used for taking samples and their preparation; manufacturer instructions for application/use must be followed.

### Diagnostics.

Diagnoses are only to be made by authorized and trained persons. Valid nomenclatures must be used. Further tests must be selected and implemented according to recognized methods.

### Storage

The reagents must be stored at +15°C to +25°C, at +5°C to +30°C. After the first opening of the bottle the contents can be used up to the expiry date when stored at +15°C to +25°C, at +5°C to +30°C.

### Shelf-life.

The bottles must be kept tightly closed at all times. The reagents must be used by the expiry date stated. Instructions for use. In order to avoid errors, the staining process must be carried out by an expert. Only for professional use. National guidelines for work safety and quality assurance must be followed. Microscopes equipped according to the standard must be used. If necessary use a standard centrifuge suitable for medical diagnostic laboratory.

### Protection against infection.

Effective measures must be taken to protect against infection in line with laboratory guidelines. Instructions for disposal. Used solutions and solutions that are past their shelf-life must be disposed of as special waste according to local disposal guidelines.

Le informazioni sopra indicate sono riportate con la massima accuratezza e rappresentano le migliori informazioni attualmente disponibili a noi. Tuttavia, non diamo garanzia di esattezza o qualsiasi altra garanzia, espressa o implicita al riguardo di tali informazioni. Inoltre; non assumiamo nessuna responsabilità derivata dal relativo uso. Gli utenti dovrebbero effettuare le loro proprie indagini per determinare l'idoneità delle informazioni per i loro scopi precisi. In nessun caso D.D.K. sarà responsabile per tutti i reclami, perdite, o danni diretti o indiretti, o verso terzi, o per i profitti persi, o danni speciali, indiretti o fortuiti, conseguenti o esemplari che possono intervenire, anche se D.D.K. si è raccomandata della possibilità di tali danni.

The information above is believed to be accurate and represents the best information currently available to us. However, we make no warranty of merchantability or any other warranty, express or implied, with respect to such information, and we assume no liability resulting from its use. Users should make their own investigations to determine the suitability of the information for their particular purposes. In no event shall D.D.K. be liable for any claims, losses, or damages of any third party or for lost profits or any special, indirect, incidental, consequential or exemplary damages, howsoever arising, even if D.D.K. has been advised of the possibility of such damages.