Wright stain stock solution technical information  
 Technical card code 09-315  
 Product code 09-315  
 Stability of product properly conserved 24 month.  
 Pack 100-200 ml or on request

Produce in Italy by:  
 DDKItalia S.r.l  
 Via Marche 19-27029 Vigevano (Pv)  
[info@ddkitalia.com](mailto:info@ddkitalia.com) • [info@dkitalia.eu](mailto:info@dkitalia.eu)

in case of emergency UE number		112
in case of emergency UK number		999
en cas d'urgence Suisse number		145

## Staining of blood and bone marrow smears and clinical-cytological specimens

### Principle

The typical color of cell nuclei, namely purple, is due to molecular interaction between eosin Y and an azure B-DNA complex. Both dyes build up the complex later. The intensity of the staining depends on the azure B content and on the ratio azure B eosin Y. The staining result can be influenced by several factors such as the pH of the solutions and buffer solution, buffer substances, fixation, staining time.

### Sample material

Air-dried blood and bone marrow smears. Clinical specimens such as urine sediment, sputum, FNAB, imprints, lavages.

### Procedure

Air-dried smears  
 Staining rack  
 Wright's solution 1 minute  
 Buffer solution (1 ml) add, mix, stain 4 minute  
 Rinse with buffer solution  
 Dry

### Staining in cuvette

Wright's solution 3 minute  
 Dilute Wright's solution 6 minute  
 Rinse with buffer solution 2 x 1 minute  
 Dry

### Result

Nuclei	red to violet
Lymphocytes	plasma blue
Monocytes	plasma gray-blue
Neutrophilic granulocytes	granules light violet
Eosinophilic granulocytes	granules brick-red to red-brown
Basophilic granulocytes	granules dark violet to black
Thrombocytes	violet
Erythrocytes	reddish

Wright stain stock solution technical information  
Technical card code 09-315  
Product code 09-315

#### Endnotes

- 1 The timing suggested in the leaflet are approximate and may vary according to your specific needs. If they are used intensively, for staining solutions may lose their dyes, so it is necessary to extend the time of staining solutions, or replace with new products.
2. Include positive control slides in each session.
3. Some hydraulic systems deliver acidic water, unsuitable for use for the part of the procedure for the blue coloration. If tap water is acidic, instead using a dilute alkaline solution, for example, water buffered by Scott.
4. The presence of purple or red-brown nuclei a blue color indicates unsatisfactory.
5. If you over-eosin staining, nuclear staining may be masked. If done correctly, with eosin staining shows a three-tone effect. To increase the differentiation of eosin, extend the time of immersion in alcohol, or use a first alcohol with a higher water content. You can adjust the times of immersion in alcohol to obtain an adequate eosin staining.
6. We do not recommend the addition of stock solution in the working solutions of haematoxylin and eosin.
7. Avoid excessive drag (carryover) of water solutions in alcoholic eosin.
8. The data generated by this procedure are to be used only to support the diagnosis and should be evaluated in conjunction with other tests and diagnostic data

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The information above is believed to be accurate and represents the best information currently available to us. However, we make no warranty of merchantability or any other warranty, express or implied, with respect to such information, and we assume no liability resulting from its use. Users should make their own investigations to determine the suitability of the information for their particular purposes. In no event shall D.D.K. be liable for any claims, losses, or damages of any third party or for lost profits or any special, indirect, incidental, consequential or exemplary damages, howsoever arising, even if D.D.K. has been advised of the possibility of such damages.