




IVD dispositivo medico-diagnostico in vitro

Shorr for sperm morphology technical information
 Technical card code 13-104
 Product code 13-104
 Pack 2x100 ml or on request
 Stability of product properly conserved at 15-20°C 24 months

Made in Italy by
 DDKItalia S.r.l
 Via Marche, 19 • 27029 Vigevano (I)
info@ddkitalia.com • www.ddkitalia.com

| | | |
|--------------------------------|---|-----|
| in case of emergency UE number |  | 112 |
| in case of emergency UK number |  | 999 |
| en cas d'urgence Suisse |  | 145 |

Preparation of semen smears

Rapid addition of fixative to semen does not permit adequate visualization of spermatozoa, as they are obscured by denatured seminal proteins. For morphological analysis, it is customary to prepare semen smears that are air-dried before fixation and staining. However, such a process leads to morphological artefacts, since air drying of semen smears is associated with: changes in sperm dimensions: dried, fixed and stained spermatozoa are smaller than live spermatozoa visualized in semen (Katz et al., 1986); expansion of immature sperm heads (Soler et al., 2000); and loss of osmotically sensitive cytoplasmic droplets, although large amounts of excess residual cytoplasm are retained. Two or more smears should be made from the fresh semen sample in case there are problems with staining or one slide is broken. Assessment is performed in replicate, preferably on each of the two slides, because there may be significant between-slide variation in sperm morphology.

Method

Mix the semen sample well

Before removing an aliquot of semen for assessment, mix the sample well in the original container, but not so vigorously that air bubbles are created.


This can be achieved by aspirating the sample 10 times into a wide-bore (approximately 1.5 mm diameter) disposable plastic pipette (sterile when necessary).

Do not mix with a vortex mixer at high speed as this will damage spermatozoa.

Remove an aliquot immediately, allowing no time for the spermatozoa to settle out of suspension. Remix the semen sample before removing replicate aliquots.

Semen smearing methods for sperm morphology

(a) "Feathering" method for undiluted semen. The semen drop (S) spreads along the back edge of the angled slide and is pulled forwards over the slide to form the smear. (b) Pipette method for washed samples. A drop of the sperm suspension (SS) is spread over the surface of the slide by pushing the horizontal pipette (P).



Normal semen samples

Clean both surfaces of the frosted slides by rubbing vigorously with lint-free tissue paper.

Label the frosted portion with identifying information (e.g. identification number, date) using a pencil. Apply a 5-10 µl aliquot of semen, depending on sperm concentration, to the end of the slide. Use a second slide to pull the drop of semen along the surface of the slide. If the dragging slide is non-frosted, the edges of both ends of the slide can be used to make four different smears. Allow the slides to dry in air and stain

IVD dispositivo medico-diagnostico in vitro

Shorr for sperm morphology technical information
 Technical card code 13-104
 Product code 13-104

Method

Fixing the air-dried semen smear

Immerse slides in acetic ethanol or 75% ethanol for 1 hour.

Sequentially immerse the slides in:

| | |
|------------------------|--------------|
| 1) Running tap water | 12-15 dips. |
| 2) Harris haematoxylin | 1-2 minutes. |
| 3) Running tap water | 12-15 dips. |
| 4) Ammoniacal ethanol | 10 dips. |
| 5) Running tap water | 12-15 dips. |
| 6) Ethanol 50% | 5 minutes. |
| 7) Shorr stain | 3-5 minutes. |
| 8) Ethanol 50% | 5 minutes. |
| 9) Ethanol 75% | 5 minutes. |
| 10) Ethanol 95% (| 5 minutes. |

*One dip corresponds to an immersion of about 1 second.

Note: The slides can be viewed unmounted or mounted.

Result

Red-stained sperm flagella are stained in samples of good and bad motility.

Both red- and blue-tailed spermatozoa are stained.

Shorr-stained slides without a cover slip revealed the presence of red and blue-stained sperm tails, but after cover slipping, significantly fewer blue-stained sperm tails can be observed.

Reagent

| | |
|---------------------|--------|
| Harris haematoxylin | 100 ml |
| Shorr | 100 ml |

References

WHO laboratory manual for the examination and processing of human semen, 5th edition, WHO, 2010

Le informazioni sopra indicate sono riportate con la massima accuratezza e rappresentano le migliori informazioni attualmente disponibili a noi. Tuttavia, non diamo garanzia di esattezza o qualsiasi altra garanzia, espressa o implicita al riguardo di tali informazioni. Inoltre; non assumiamo nessuna responsabilità derivata dal relativo uso. Gli utenti dovrebbero effettuare le loro proprie indagini per determinare l'idoneità delle informazioni per i loro scopi precisi. In nessun caso D.D.K. sarà responsabile per tutti i reclami, perdite, o danni diretti o indiretti, o verso terzi, o per i profitti persi, o danni speciali, indiretti o fortuiti, conseguenti o esemplari che possono intervenire, anche se D.D.K. si è raccomandata della possibilità di tali danni.

The information above is believed to be accurate and represents the best information currently available to us. However, we make no warranty of merchantability or any other warranty, express or implied, with respect to such information, and we assume no liability resulting from its use. Users should make their own investigations to determine the suitability of the information for their particular purposes. In no event shall D.D.K. be liable for any claims, losses, or damages of any third party or for lost profits or any special, indirect, incidental, consequential or exemplary damages, howsoever arising, even if D.D.K. has been advised of the possibility of such damages.