

IVD dispositivo medico-diagnostico in vitro

Leuco test technical's information's

Technical's card 13-106

Product code 13-106

Pack 2x50 ml or on request

Stability of product properly conserved at 4°C 12 months

CND code W01030799

Producers:

DDKItalia S.r.l

Via Marche, 19 • 27029 Vigevano (Pv)

info@ddkitalia.com•www.ddkitalia.com

For emergency contact your nearest anti-poison unit.

Principle

Most human ejaculates contain leucocytes, the predominant cell type being the neutrophil. Excessive presence of these cells (leucocytospermia) may indicate the existence of reproductive tract infection.

Furthermore, leucocytospermia may be associated with defects in the semen profile, including reductions in the volume of the ejaculate, sperm concentration, and sperm motility, as well as loss of sperm function as a result of oxidative stress (Aitken et al, 1989; Aitken and West, 1990) and - or secretion of cytotoxic cytokines (Hill et al, 1987).

It is difficult to define a threshold concentration of leucocytes beyond which fertility will be impaired. The impact of these cells depends upon the site at which the leucocytes enter the semen, the type leukocyte involved, and their state of activation. As a general rule, a normal ejaculate should not contain more than 5×10^6 round cells/ml, while the number of leucocytes should not exceed 1×10^6 /ml (WHO, 1992). When the semen contains more than 1×10^6 /ml white blood cells, microbiological tests should be performed to investigate if there is an accessory gland infection.

The granules in the neutrophilic polymorphous leucocytes contain peroxidase, which together with the hydrogen peroxide form water and free oxygen ions, these oxygen ions oxidize the benzidine which colours brown, thus staining the cells brown. Eosin 10b 1 contrast fluid to differentiate peroxidase positive round cells from peroxidase negative round cells.

Note. The absence of leucocytes does not exclude the possibility of an accessory gland infection.

Method

1) Preparation of work solution:

Add 30 μ l of Hydrogen peroxide (reagent B) to 1 ml of benzidine (reagent A). This working solution remains stable for 1 day.

2) Mix 1 drop (10 μ l) of sperm with 1 drop (10 μ l) of working solution, using the edge of the cover slip.

Mix thoroughly for at least 1 minute.

3) Cover with the cover slip approximately two minutes after initial mixing, avoiding air bubbles.

Formation of small air-bubbles is normal and due to peroxidase reaction. The higher the concentration of peroxidase positive cells, the more bubbles will form.

4) Read the result after 2 minutes at a magnification of 400x (read at least 20 separate microscope fields).

Calculation of the concentration of white blood cells

Known concentration of spermatozoa. Count the number of WBC and the number of spermatozoa. Calculate the concentration of WBC based on this formula. (Number of WBC - number of sperm cells) x sperm concentration (mill/ml). This method will only work if the semen sample contains sperm cells (preferably more than 10 mill/ml).

Unknown concentration of spermatozoa

In this case the concentration of WBC can be determined by multiplying the number of WBC with a known factor based on the size of a microscope field and the height between the object glass and the cover glass (or the depth of the semen sample).

The diameter of the microscope field can be measured using a micrometer. The surface area in one field is equal to the square of the radius multiplied with pi ($s = \pi r^2$).

e.g. Diameter = 250 μ m \Rightarrow radius = 125 μ m \Rightarrow surface area = 49086 μ m²

IVD dispositivo medico-diagnostico in vitro

Leuco test technical's information's
Technical's card 13-106
Product code 13-106

The distance between the object glass and the cover glass can be calculated using the formula:
height (in μm) = volume (in μl) / (length x width of cover glass in mm) e.g. volume of the sample = 20 μl ,
coverglass = 24 x 40 mm
Height = 20 / (24 x 40) = 0.0208mm or 20.8 μm
Knowing these figures a factor can be determined using this formula:
Factor = 1 000 000 / (surface area x height)
e.g. Factor = 1 000 000 μm^3 / (49086 μm^2 x 20.8 μm) = 0.98
This means that if 5 WBC's are counted in one microscope field, the corresponding concentration would be 4.9 mill/ml.

Result

Yellow to brown stained cells are peroxidase positive cells: neutrophil polymorphous leucocytes
Pink stained cells: all other cells

Limitations of the method

LeucoScreen only stains peroxidase positive white blood cells, other types of white blood cells (e.g. lymphocytes and monocytes) cannot be detected.

Reagent

Benzidine (WHO modified) solution A	50 ml
Hydrogen peroxide solution B	50 ml

Remarks

Formation of a sediment in Reagent 1 is normal. Simply pour Reagent 1 over filter-paper to eliminate sediment.

References

WHO laboratory manual for the examination and processing of human semen, 5th edition, WHO, 2010

Le informazioni sopra indicate sono riportate con la massima accuratezza e rappresentano le migliori informazioni attualmente disponibili a noi. Tuttavia, non diamo garanzia di esattezza o qualsiasi altra garanzia, espressa o implicita al riguardo di tali informazioni. Inoltre; non assumiamo nessuna responsabilità derivata dal relativo uso. Gli utenti dovrebbero effettuare le loro proprie indagini per determinare l'idoneità delle informazioni per i loro scopi precisi. In nessun caso D.D.K. sarà responsabile per tutti i reclami, perdite, o danni diretti o indiretti, o verso terzi, o per i profitti persi, o danni speciali, indiretti o fortuiti, conseguenti o esemplari che possono intervenire, anche se D.D.K. si è raccomandata della possibilità di tali danni.

The information above is believed to be accurate and represents the best information currently available to us. However, we make no warranty of merchantability or any other warranty, express or implied, with respect to such information, and we assume no liability resulting from its use. Users should make their own investigations to determine the suitability of the information for their particular purposes. In no event shall D.D.K. be liable for any claims, losses, or damages of any third party or for lost profits or any special, indirect, incidental, consequential or exemplary damages, howsoever arising, even if D.D.K. has been advised of the possibility of such damages.