

Sample kit MCDh

REF. 960200-0000

Fixation and differential staining of cellular structures



IFU007A-RAL

For professional use only.

Please read all this information carefully before using this device.

Table of contents

Intended use.....	1
Principle.....	1
Kit description	2
Storage	2
Active components	2
Hazard classification and safety information.....	3
Personnel qualification	3
Specific equipment and reagents required but not provided	3
Operating procedure.....	4
Expected results.....	5
Performance	6
User quality control	6
Other products.....	6
Recommendations, notes, and troubleshooting.....	6
Table of symbols and abbreviations	8
Bibliography.....	8
Changes tracking.....	8

Intended use

MCDh staining is intended to be used for the fixation and the differential staining of cellular structures prior microscopic examination.

If applicable, RAL Diagnostics recommends using the associated RAL Diagnostics products and cannot guarantee that the expected results will be achieved if used in combination with products of other brands.

Principle

The MCDh panoptic staining allows to perform blood cell counting, realized by using successively four reagents: MCDh 1, MCDh 2, MCDh 3 and MCDh 4.

As a variation of May-Grünwal Giemsa staining, MCDh staining also enables a cytological study of fluids (urines, cerebrospinal, articular, pleural, etc.) and solid samples (biopsies, tissues, swabs etc.). The interpretation of smears provides qualitative and quantitative evaluation of the presence of bacteria, leukocytes, and epithelial cells.

MCDh 1, formulated with ethylic alcohol, is a mixture of neutral stains. It allows a smear fixation and prepares the staining, especially the one of hydrosoluble elements such as basophilic granules. Those stains are inactive in alcoholic medium, and only react selectively when released in aqueous MCDh 2 solution. This releasing generates the precipitation of neutral stains, leading to erythrocytes, cytoplasm of some leukocytes like neutrophilic granulocytes, cytoplasm of bacteria as well as eosinophilic granules staining. MCDh 3 is an aqueous solution which stains cytoplasm of monocytes and lymphocytes. MCDh 3 also eases the metachromasia process as it colors azurophilic granules red.

Eventually, MCDh 4 removes excess of stain and participates to differentiation of cellular elements thanks to action of specially selected rinsing agents.

The successive action of MCDh 1, MCDh 2, MCDh 3 and MCDh 4 brings the violet color (typical Romanowsky-Giemsa effect), particularly visible in chromatin, platelets, and neutrophilic granules.

Kit description

MCDh 1

Clear dark blue solution
REF. 313590-0100 1 X 100 mL

MCDh 2

Clear colorless solution
REF. 313570-0100 2 X 100 mL

MCDh 3

Clear dark blue solution
REF. 313560-0100 1 X 100 mL

MCDh 4


Clear colorless solution
REF. 313600-0100 1 X 100 mL

For a specific batch, refer to the analysis certificate of the batch available at my.ral-diagnostics.fr.

Storage

Storage temperature: 15-25°C away from light.

Bottle shelf life before and after opening: refer to the expiry date on the label.



Active components

MCDh 1

May-Grünwald: ca 0.1%
Methylene azure I blue – CAS - 531-55-5: ca 0.05%

MCDh 2

Potassic mono phosphate - CAS 7778-77-0: ca 0.05%
Anhydrous disodic phosphate - CAS 7558-79-4: ca 0.04%

MCDh 3

Methylene blue – CAS - 61-73-4: < 0.25%

MCDh 4

Potassic mono phosphate - CAS 7778-77-0: ca 0.03%
Anhydrous disodic phosphate - CAS 7558-79-4: ca 0.03%

Hazard classification and safety information

MCDh 1

Danger: H225-Highly flammable liquid and vapour.

P210-Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.



MCDh 2

Warning: H317-May cause an allergic skin reaction. H412 - Harmful to aquatic life with long lasting effects.

P261 - Avoid breathing dust/fume/gas/mist/vapours/spray. P280 - Wear protective clothing, protective gloves, eye protection. P333+P313 - If skin irritation or rash occurs: Get medical advice/attention. P362+P364 - Take off contaminated clothing and wash it before reuse.



CONT	5-chloro-2-methyl-2H-isothiazol-3-one 2-methyl-2H-isothiazol-3-one
-------------	---

MCDh 3

No labelling applicable

MCDh 4

Warning: H226 - Flammable liquid and vapour. H317 - May cause an allergic skin reaction. H412 - Harmful to aquatic life with long lasting effects.

P210 - Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P261 - Avoid breathing dust/fume/gas/mist/vapours/spray. P280 - Wear protective gloves, protective clothing, eye protection. P333+P313 - If skin irritation or rash occurs: Get medical advice/attention. P362+P364 - Take off contaminated clothing and wash it before reuse.



CONT	5-chloro-2-methyl-2H-isothiazol-3-one 2-methyl-2H-isothiazol-3-one
-------------	---

Personnel qualification

All samples and products must be handled by qualified and authorized personnel, using individual or collective protection, in accordance with the national directives in force in the laboratories. Personnel must also be aware of the classification of hazardous materials indicated on the label and the safety data sheet (available at my.ral-diagnostics.fr).

The specimen must be treated in accordance with procedures available in the laboratory and required by national authorities.

The diagnosis must be conducted by qualified and authorized personnel, in accordance with the procedures in force within the laboratory.

Specific equipment and reagents required but not provided

Microscope slides and absolute ethanol.

This equipment may vary depending on the protocol. Please refer to the relevant protocol (see the section operating procedure) to ensure that you have the necessary equipment to carry out tests.

Operating procedure

The equipment used for sample processing must comply with the supplier's instructions for use.

Sample preparation

Manual blood smear: Mix the tube by slow inversion and install a smearing droplet device. Invert the tube and lightly press the drop depositor onto a slide to deposit a small drop of blood (Fig. 1- slide A at step 1).

Using another slide tilted at 45° (Fig. 1- slide B at step 1), spread the blood by capillarity on the short edge (Fig. 1- steps 2 & 3) using a pushing motion (Fig. 1- step 4). A good quality smear does not reach the end of the slide and has a gradual decrease in thickness until the end is feathered. Allow the smear to air dry before fixing or staining.

NB: if you do not have a smearing droplet device, open the tube, and use a pipette to deposit a blood drop.




Figure 1. Schematic representation of performing a blood smear

A & B: slides, 1 – 4: steps 1 to 4

Manual bacterial smear: made a thin film of bacteria sample and leave the slide to dry at room temperature. Then bacterial smear can be fixed with mild heat

source (Bunsen burner or hot plate) or chemically fixed with chemical fixative (methanol, ethanol, acetic acid, or formalin...)

NB: Never pass through the flame a smear that is not entirely dry, this could cause cracklings and dissemination of bacteria (creation of aerosols).

If necessary, the two fixations can be combined.

Reagents and instruments preparation

No preparation needed. The solutions are ready to use and the reagents containers have been designed to be used for slides staining.

Protocols

The staining steps of the protocols indicated below consist of a successive dipping of the slides in the different staining baths.

Protocol for hematological samples - Manual bath staining method - CellaVision® automates analysis

Processing time: 11 min 30

Steps	Reagent	Time [mm: ss]	Indications
Fix and pre-stain	MCDh1	06: 00	Do not agitate
Stain	MCDh2	01: 00	Agitate in the bath 5 to 10 times at the end of countdown
Stain	MCDh2	02: 00	Agitate in the bath 5 to 10 times at the beginning and the end countdown
Stain	MCDh3	00: 30	
Rinse	MCDh4	02: 00	
Dry	No	≥03: 00	No

Protocol for hematological samples- Manual bath staining method - Manual microscopic analysis

Processing time: 09 min 50

Steps	Reagent	Time [mm: ss]	Indications
Fix and pre-stain	MCDh1	06: 00	Do not agitate
Stain	MCDh2	01: 00	Agitate slowly 3 times in the bath at the end of countdown
Stain	MCDh2	02: 00	
Stain	MCDh3	00: 40	
Rinse	MCDh4	00: 10	Agitate continuously in the bath during countdown
Dry	No	≥03: 00	No

Note: In case of refringence/water artefact phenomena, pre-fix the slides 2min in a bath of absolute ethanol before staining. Directly start the staining after pre-fixation step without drying the slides.

Protocol for cyto-bacterial samples- Manual bath staining method- Manual microscopic analysis

Processing time: 09 min 32

Steps	Reagent	Time [mm: ss]	Indications
Fix and pre-stain	MCDh1	07: 00	Do not agitate
Stain	MCDh2	01: 00	Agitate slowly in the bath at the end of countdown
Stain	MCDh2	01: 00	
Stain	MCDh3	00: 30	
Rinse	MCDh4	00: 02	Agitate continuously in the bath during countdown
Dry	No	≥03: 00	No

Expected results

Hematological samples:

- Nuclei / chromatin:** +/- dense purple
- Granulocytes cytoplasm:** light purplish-pink
- Granulocytes eosinophilic granules:** orangey
- Granulocytes basophilic granules:** dark blue
- Granulocytes neutrophilic granules:** +/- deep purple
- Lymphocytes cytoplasm with RNA:** pure blue
- Lymphocytes cytoplasm without RNA:** light blue
- Lymphocytes azurophilic granules:** red
- Monocytes cytoplasm:** cloudy blue
- Erythrocytes:** pinkish-beige
- Platelets chromomere:** purplish-red
- Platelets hyalomere:** bluish
- Blood parasites nucleus:** red
- Blood parasites cytoplasm:** blue

Cyto-bacterial samples

- Nuclei / chromatin:** +/- dense purple
- Cytoplasm:** light blue to light purplish- blue
- Erythrocytes:** pinkish-beige
- Bacteria:** blue to purplish-pink

If observed results vary from those expected, please contact RAL Diagnostics technical service through your usual supplier for assistance.

Performance

This medical device is state of the art. Its analytical performance, scientific validity and medical relevance are assessed in the CE marking review.

To ensure product performance, use clean and dry laboratory equipment.

The laboratory is responsible for notifying the manufacturer and state competent authority of any serious incident relating to the medical device uses.

User quality control

Users are responsible for determining the appropriate quality control procedures for their laboratory and complying with applicable laboratory regulations.

Hematological sample: RAL Diagnostics recommends staining a freshly made blood smear with a normal WBC count and no known abnormal pathology at reagent renewal and for the first staining cycle each day. Slides stained for quality control purposes should be checked to ensure that they are satisfactory for intended test (properly stained and free of precipitate).

Bacterial sample: RAL Diagnostics recommend using a known bacteria sample for reagents quality control at reagents renewal, for each staining cycle or at least for the first staining cycle if a stain is performed multiple times daily. The result is checked under a microscope, in comparison with the results obtained by the usual technique validated by the laboratory. Staining results must also be compliant with this manual expected results.

These quality control procedures should only be performed by qualified personnel.

Other products

For more information contact your usual supplier.

Recommendations, notes, and troubleshooting

Products appearance

If the appearance of the products differs from that describe above in this manual, do not use it and contact RAL Diagnostics technical service for assistance.

Procedures notes

To prevent any degradation of the products, please respect the storage and handling conditions specified in this manual.

In case of refringence/water artefact phenomena, pre-fix the slides 2min in a bath of absolute ethanol before staining.

Products stability

Each RAL Diagnostics product can be used until the expiration date indicated on, in its original packaging and always hermetically sealed.

Staining stability

Staining quality and reproducibility depend on the correct use of the products. Staining conducted according to these recommendations will remain stable for several days. If it is necessary to store the stained smear for several months or years, RAL Diagnostics recommends mounting them with a coverslip using a suitable mounting liquid and storing them in a light and dustproof container.

Instructions for cleaning and waste disposal

All biological samples, effluents and used consumables should be considered potentially hazardous.



To avoid any risk, apply the following instructions: dispose of samples, effluents and consumables in accordance with laboratory standards and applicable national and local standards and regulations.

Chemical and biological waste must be collected and processed by specialized, registered companies.

Table of symbols and abbreviations

Depending on the product, you may find the following symbols on the device or the packaging material.

GHS PICTOGRAMS	INTERPRETATION	SYMBOL	INTERPRETATION
	Explosive		Batch code
	Flammable		Serial number
	Oxidizer		Catalogue reference
	Compressed gas		Date of manufacture
	Corrosive		Use up to
	Toxic		Unique device identifier
	Harmful		Manufacturer
	Health Hazard		Importer
	Environmental Hazard		Entity distributing the medical advice in the region concerned
	No labelling applicable		CE marking device
			In vitro diagnostic medical device
			Authorized Representative in the European Community
			Authorized Representative in Switzerland
			Complies with UK guidelines
			UK CA
			Do not use if packaging is damaged
			Keep away from light
			Temperature limit: 15-25°C
			Temperature limit: 15-30°C
			Keep dry
			Box: handling upwards
			Fragile
			Sterilised by irradiation
			Single sterile barrier system with outer protective packaging
			Sterile and radiation-sterilised barrier suit
			Do not reuse
			Do not sterilize
			Contents sufficient for n tests
			Hazardous material contained
			Consult instructions for use
			Use
			After opening, use within XX months
			The product must not be used in conjunction with an automatic colouring machine
			Indicates a medical device that contains potentially carcinogenic, mutagenic or reprotoxic (CMR) substances, or substances classified as endocrine disruptors

Bibliography

BENATTAR L., FLANDRIN G., *Morphometry and Quality Control for a May-Grünwald Giemsa stained preparation. A 40 centers cooperative study. Leuk. & Lymphoma* 1999, 33, 587-591.

BENATTAR L., FLANDRIN G., *Etapes de l'automatisation de l'étude microscopique du sang. Rencontre Médecins biologistes, 2002. ATEB, Journée Technique Parisienne, mars 1977.*

DUHAMEL G., DUHAMEL E., *Cytologie hématologique, Les cellules pathologiques I et II, Coloration au May-Grünwald Giemsa RAL, Biologiste et Praticien et Réactifs RAL, 1984 et 1989.*

Ecole Nationale de Chimie, *Coloration de Pappenheim, Présentation théorique des mécanismes cytochimiques des colorants neutres avec applications techniques détaillées, Journée du technicien biologiste, mars 1980, p. 1-9.*

GENTILHOMME O., TREILLE-RITOUET D., BRYON P-A., *Cytologie hématologique, Les cellules normales, Coloration au May-Grünwald Giemsa RAL, Réactifs R.A.L, 1989.*

THEML H., *ATLAS de poche d'Hématologie, Médecine-Sciences Flammarion, p. 19-25, 2000*

Changes tracking

Date	Version	Changes
05/2022	IFU007A-RAL	IVDR (EU) 2017/746 compliance