

Schiff reagent

REF.320680

Multi-application microscopy reagent



IFU123A

Changes tracking 7

Legal representatives 7

For professional use only.

Please read all information carefully before using this device.

IFU content may change, make sure you have the latest version available at my.ral-diagnostics.fr.

Table of contents

Intended use.....	1
Principle.....	1
Device description	2
Storage and use conditions	2
Hazard classification and safety information.....	2
Personnel qualification	2
Specific equipment and reagents required but not provided	2
Operating procedure.....	3
Expected results.....	5
Performance	6
User quality control	6
Other products.....	6
Recommendations, notes and troubleshooting.....	6
Table of symbols and abbreviations	7
Bibliography.....	7

Intended use

Schiff reagent is intended to be used as multi-application microscopy reagent particularly in histo-cytological application prior microscopic examination.

If applicable, CellaVision RAL Diagnostics recommends using the associated CellaVision RAL Diagnostics products and cannot guarantee that the expected results will be achieved if used in combination with products of other brands.

Principle

Schiff reagent in combination with other staining devices allows histo-cytology and mycology staining including Periodic Acid Schiff (PAS) and Feulgen staining. PAS stains carbohydrates and glycoproteins in the wall of some pathogenic factors and helps to detect some parasitic and fungal elements. Its combine periodic acid as oxidizing agent, that breaks the bonds between two carbons of some chemical groups and Schiff reagent that stain released aldehydes functions. PAS can be applied on smears (cutaneous squama, oropharynx and esophagus cells, broncho-alveolar lavage fluid, smears through affixing of organs), or on histological sections.

Schiff reagent is also used for DNA assay in histo-cytological samples by two successive steps partial hydrolysis by HCl removing the purine bases (adenine and guanine), then this apuric DNA is stained by Schiff reagent.

Device description

Schiff reagent

Clear colorless solution

REF. 320680-0500 1 X 0.5 L

REF. 320680-1000 1 X 1 L

For a specific batch, refer to the analysis certificate of the batch available at my.ral-diagnostics.fr.


Storage and use conditions

Storage and use temperature: 15-25°C.

Storage and use conditions: away from light and heat sources.

Bottle shelf life before opening: refer to expiry date on the label.

Bottle shelf life after opening : refer to expiry date on the label and if the "period after opening" symbol is present take it into account.



Hazard classification and safety information

Schiff reagent

Warning:

H319 - Causes serious eye irritation.

P280 - Wear protective gloves, protective clothing, eye protection.

P337+P313 - If eye irritation persists: Get medical advice/attention.



Personnel qualification

All samples and products must be handled by qualified and authorized personnel, using individual or collective protection, in accordance with the national directives in force in the laboratories. Personnel must also be aware of the classification of hazardous materials indicated on the label and the safety data sheet (available at my.ral-diagnostics.fr).

The diagnosis must be conducted by qualified and authorized personnel, in accordance with the procedures in force within the laboratory.

Specific equipment and reagents required but not provided

Hydrochloric acid, potassium metabisulphite, periodic acid, acetic acid, absolute ethanol, acetone, ether, sodium metabisulfite, lithium carbonate, microscope slides and these following RAL Diagnostic reagents:

Mayer Haemalum REF. 320550

Saffron in alcoholic solution REF. 369200

This equipment may vary depending on the protocol. Please refer to the relevant protocol (see the section operating procedure) to ensure that you have the necessary equipment to carry out tests.

Operating procedure

The equipment used for sample processing must comply with the supplier's instructions for use.

Sample preparation

The specimen must be treated in accordance with procedures available in the laboratory and required by national authorities.

Histological sections: dewax and et hydrate tissues sections in appropriate reagents before staining.

Reagents and instruments preparation

The Schiff reagent is ready to use.

Sulphureous acid: Prepare extemporaneously 5 ml of HCl N + 5 ml of 10%-Potassium Metabisulphite + 100 ml of distilled water.

Periodic acid solution: prepare between 0.5 and 0.8% periodic acid solution.

Protocols

The staining steps of the protocols indicated below consist of a successive covering of the slides with the different staining reagents or dipping of the slides in the different staining baths. Please refer to the title to know which case you are in. For the covering method, place slide on a stand with fixed smear on top. The processing time only considers the dipping time in the reagents.

Protocol for Feulgen Staining- Bath staining method - Manual microscopic analysis

Processing time: 01:46:00

Steps	Reagent	Time [mm: ss]	Indications
Rinse	HCl 1N solution	01 :00	Optional
Hydrolyze	HCl 1N solution	05 :00	Can be extended to 20 min according to fixative*
Stop hydrolyzing	Cold distilled water	NA	Stop hydrolyzing rinsing brusquely with cold distilled water
Stain	Schiff reagent	90:00	Can be extended to 2 hours then dry
Discolor	Sulphureous acid	15:00	Dip in 3 successive baths of 5 minutes in each
Rinse	Distilled water	NA	NA
Dehydrate	Croissant degrees ethanol	NA	To absolute ethanol
Dehydrate	Xylene or toluene	NA	Pass slide in
Mount	Xylene or toluene mounting media	NA	NA

Note: Tissues sections are preferably fixed in Zenker, Helly, Baker, Carnoy, sublimate formalin or sublimate ethanol.

* See underlying table.

Specificity Control: Carry out the same reaction on another section cancelling hydrolysis (steps 1 to 3, Rinse, hydrolyze and stop hydrolyzing): DNA should not be colored.

Fixative	Time in min	Fixative	Time in min
Bouin-Allen	22	Helly	8
Sublimate Bouin-Allen	14	Stowell	20
Carnoy	8	Flemming	16
Formalin	8	Champy	25
Sublimate Formalin	8	Susa	18
Absolute ethanol	5	Regaud	14
Acetic Formalin Alcohol	7	Newcomer	20
Zenker	5		

Table of hydrolysis time at 60°C in HCl N according to the fixative

Protocol for periodic acid Schiff staining in Mycology and Parasitology - Manual covering method - Manual microscopic analysis

Processing time: 02:00:30

Steps	Reagent	Time [mm: ss]	Indications
Fix	30% acetic acid aqueous solution	NA	Dry heating it very gently until total evaporation of the solution
Fix	50/50 mixture of ethanol and acetone	NA	Dry until total evaporation of the solution, repeat 3 times
Fix	Ether	NA	Dry until total evaporation of the solution, repeat 2 times
Fix	70° ethanol	15:00	Can be extended to 30 min
Fix	95° ethanol	15:00	Pass in a bath
Rinse	Running water	NA	Running water
Pre-stain	Periodic acid	10:00	NA
Rinse	Running water	NA	NA

Stain	Schiff reagent	16:00	NA
Develop stain	Sodium metabisulfite	06:00	Without rinsing previous reagent
Rinse	Running water	05:00	NA
Stain	Mayer haemalum	02:30	NA
Rinse	Running water	No	Until turn in color
Differentiate	Lithium carbonate	15:00	saturated and filtered aqueous solution
Rinse	Running water	NA	NA
Rinse	Absolute ethanol	NA	NA
Stain	Saffron in alcoholic solution	01: 00	NA
Rinse	Absolute ethanol	20: 00	2 x 10 min Can be extend to 30 min (3 x 10 min)
Dehydrate	Xylene or toluene	15: 00	NA
Mount	Xylene or toluene mounting media	NA	NA

Protocol for periodic acid Schiff staining of tissues sections in Histology - Manual bath method - Manual microscopic analysis

Processing time: 00:39:00

Steps	Reagent	Time [mm: ss]	Indications
Pre-stain	Periodic acid	05:00	Precise time requirement
Rinse	Running water	05:00	Can be extended to 10 min
Rinse	Distilled water	02:00	Can be extended to 3 min
Stain	Schiff reagent	15:00	Can be extended to 30 min
Rinse	Distilled water	02:00	Can be extended to 3 min
Rinse	Running water	02:00	Can be extended to 5 min
Stain	Mayer haemalum	05:00	NA
Rinse	Running water	03: 00	Rinse and let in a bath of running water for 3 to 5 min
Dehydrate	Croissant degrees ethanol	NA	To absolute ethanol
Dehydrate	Xylene or toluene	NA	Pass slide in
Mount	Xylene or toluene mounting media	NA	NA

Expected results

Feulgen staining

DNA: pink-red

P.A.S. staining in Mycology and Parasitology

P.A.S. reaction stains red more or less bright:

- the wall of fungal mycelian and levuriform elements
- the precystic and cystic forms of *Pneumocystis Carinii*
- some protozoa: trophozoites of *Entamoeba histolytica*, as well as the forepart of the *microsporidia* spore filament

Concerning metazoa, P.A.S. reaction eases the detection and the characterization of the cuticular (or lamellar) membrane of vesicular larva of *Ecchinococcus granulosus* (hydatidosis) and *Ecchinococcus multicularis* (alveolar echinococcosis).

P.A.S. reaction is also very useful to the identification of the cuticle of tissular nematoda (e.g., *Anisakis* larva) and of arthropides integumentary structures (*Demodex*, *Sarcoptes*, cutaneous mites causing myiasis; *Tunga penetrans* contributing to sarcopsyllosis or similar organisms, pentastomes).

P.A.S. staining (Periodic Acid Schiff) of tissue sections in Histology

Nuclei: blue

Glucidic substances: bright red

Main substances colored by Periodic Acid Schiff are:

- Polysaccharides: Glycogen, starch, cellulose, dextrans.
- Glycoproteins: Fucomucins of mucus cells, glycoproteins of basal cells, reticulin, collagen (slightly tinted), crystalline lens, bacterial capsules.
- Acid mucopolysaccharides: Connective fundamental substance (slightly colored), mucus cells (neither hyaluronic acids, nor chondrotin acids are stained, except some slightly sulphated).

If observed results vary from those expected, please contact CellaVision RAL Diagnostics technical service through your usual supplier for assistance.

Performance

This medical device is state of the art. Its analytical performance, scientific validity and medical relevance are assessed in the CE marking review.

To ensure product performance, use clean and dry laboratory equipment.

The laboratory is responsible for notifying the manufacturer and state competent authority of any serious incident relating to the medical device uses.

User quality control

Users are responsible for determining the appropriate quality control procedures for their laboratory and complying with applicable laboratory regulations.

These quality control procedures should only be performed by qualified personnel.

Other products

For more information, please contact your usual supplier.

Recommendations, notes and troubleshooting

Products appearance

If the appearance of the products differs from the description above, do not use it and contact CellaVision RAL Diagnostics technical service through your usual supplier for assistance.

Procedure notes

To prevent products degradation, please comply with the storage and handling recommendations specified in this manual.

In the case of Feulgen staining, to assert the positivity of the reaction, Schiff's reagent must stain the sample only after lysis and does not stain the same unlysed sample.

Product stability

Every CellaVision RAL Diagnostics product can be used until the expiry date indicated on, in its original packaging if it is still hermetically sealed.

Staining stability

Staining quality and reproducibility depend on the correct use of the products. Staining conducted according to these recommendations will remain stable for several days. If it is necessary to store the stained smears for several months or years, CellaVision RAL Diagnostics recommended mounting them with a coverslip, using a suitable mounting liquid and storing them in a light and dustproof container.

Instructions for cleaning and waste disposal

All biological samples, effluents and used consumables should be considered potentially hazardous.



To avoid any risk, apply the following instructions: dispose of samples, effluents and consumables in accordance with laboratory standards and applicable national and local standards and regulations.

Chemical and biological waste must be collected and processed by specialized, registered companies.

Table of symbols and abbreviations

Depending on the product, you may find the following symbols on the device or the packaging material.

GHS Pictograms	Interpretation
	Explosive
	Flammable
	Oxidizer
	Compresses gas
	Corrosive
	Toxic
	Harmful
	Health Hazard
	Environmental Hazard
	No labelling applicable

Symbols	Interpretation
	Batch code
	Serial number
	Catalogue reference
	Date of manufacture
	Use up to
	Unique device identifier
	Manufacturer
	Importer
	Entity distributing the medical advice in the region concerned
	CE marking device
	In vitro diagnostic medical device
	Authorised Representative in the European Community
	Authorised Representative in Switzerland
	Complies with UK guidelines
	Do not use if packaging is damaged
	Keep away from light
	Temperature limit: 15-25°C
	Temperature limit: 15-30°C
	Keep dry
	Box: handling upwards
	Fragile
	Sterilised by irradiation
	Single sterile barrier system with outer protective packaging
	Sterile and radiation-sterilised barrier suit
	Do not reuse
	Do not resterilize
	Contents sufficient for n tests
	Hazardous material contained
	Consult instructions for use
	Use
	After opening, use within XX months
	The product must not be used in conjunction with an automatic colouring machine
	Indicates a medical device that contains potentially carcinogenic, mutagenic or reprotoxic (CMR) substances, or substances classified as endocrine disruptors

Bibliography

GANTER P., JOLLES G., *Histochimie normale et pathologique*, ed. GAUTHIER-VILLARS, vol. 2, 1970, p. 1526-1527 and 1533-1534.

SEGRETAIN G., DROUHET E., MARIAT F., *Diagnostic de laboratoire en mycologie médicale*, Maloine, 3^{ème} éd., 1974, p. 125-126.

Changes tracking

Date	Version	Changes
02/2023	IFU123A	IVDR (EU) 2017/746 compliance

Legal representatives

Country	Address
United Kingdom	QAvis UK Ltd, company N° SC679796, 56-66 Frederick Street Edinburgh, EH21LS, United Kingdom
Switzerland (CH-REP)	MedEnvoy Switzerland, Gotthardstrasse 28, 6302 Zug Switzerland