

# Kit Trichrome de Masson, light green variation



REF. 361350-0000

IFU084A-RAL

Differential staining of histo-cytological structures

For professional use only.

Please read all information carefully before using this device.

## Table of contents

Intended use.....	1
Principle .....	1
Kit description .....	2
Storage.....	2
Hazard classification and safety information.....	3
Personnel qualification .....	3
Specific equipment and reagents required but not provided.....	3
Operating procedure .....	4
Expected results.....	5
Performance.....	5
User quality Control .....	5
Other products.....	5
Recommendations, notes, and troubleshooting.....	5
Table of symbols and abbreviations .....	7
Bibliography.....	7
Change tracking .....	7

## Intended use

Kit Trichrome de Masson, light green variation is intended to be used for differential staining of histo-cytological structures prior microscopic examination

If applicable, RAL Diagnostics recommends using the associated RAL Diagnostics products and cannot guarantee that the expected results will be achieved if used in combination with products of other brands.

## Principle

Kit Trichrome de Masson, light green variation allows a Trichromic staining of histological sections combining three stains: a nuclear stain, Mayer Haemalum, a cytoplasmic stain, a mix of acid stains: Ponceau Fuchsin and an elective stain of collagen, another acid stain: Light green.

## Kit description

### Mayer haemalum

Clear red violet solution  
REF. 3205501B100

1 X 100 mL

### Light green

Clear green solution  
REF. 361700-0100

1 X 100 mL

### Phosphomolybdic acid

Clear light yellow solution  
REF. 3617603B100

1 X 100 mL

### Ponceau fuchsin

Clear red solution  
REF. 3617702B100

1 X 100 mL

### 1% acetic water

Clear colorless solution  
REF. 3618105B100

1 X 100 mL

### 1% acetic water

Clear colorless solution  
REF. 3618406B100

1 X 100 mL

## Storage

Storage temperature: 15-25°C away from light.

Bottle shelf life before and after opening: refer to expiry date on label.



For a specific batch, refer to the analysis certificate of the batch available at [my.ral-diagnostics.fr](http://my.ral-diagnostics.fr).

## Hazard classification and safety information

### Mayer haemalum

Warning: H226 - Flammable liquid and vapour. H302 - Harmful if swallowed. H371 - May cause damage to organs. P210 - Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P264 - Wash hands thoroughly after handling. P308+P311 - IF exposed or concerned: Call a POISON CENTER or doctor.



CONT CH3OH

### Light green

No labelling applicable

### Phosphomolybdic acid

No labelling applicable

### Ponceau fuchsin

No labelling applicable

### 1% acetic water

No labelling applicable

### 1% acetic water

No labelling applicable

## Personnel qualification

All samples and products must be handled by qualified and authorized personnel, using individual or collective protection, in accordance with the national directives in force in the laboratories. Personnel must also be aware of the classification of hazardous materials indicated on the label and the safety data sheet (available at [my.ral-diagnostics.fr](http://my.ral-diagnostics.fr)).

The specimen must be treated in accordance with procedures available in the laboratory and required by national authorities.

The diagnosis must be conducted by qualified and authorized personnel, in accordance with the procedures in force within the laboratory.

## Specific equipment and reagents required but not provided

Microscope slides, alcohol baths, Toluene or Xylene or substitute, Toluene or Xylene or substitute based mounting medium.

This equipment may vary depending on the protocol. Please refer to the relevant protocol (see the section operating procedure) to ensure that you have the necessary equipment to carry out tests.

## Operating procedure

The equipment used for sample processing must comply with the supplier's instructions for use.

## Sample preparation


Dewax and et hydrate tissues sections in appropriate reagents before staining.

## Reagents and instruments preparation

No preparation needed. The solutions are ready to use and the reagents containers have been designed to be used for slides staining.

## Protocols

The staining steps of the protocols indicated below consist of a successive dipping of the slides in the different staining baths.



## Protocol for histological sections staining - Manual bath method - Manual microscopic analysis

Dewax and et hydrate tissues sections in appropriate reagents before staining.

Processing time: 28 min

Steps	Reagent	Time [mm:ss]	Indications
Stain	Mayer haemalum	10:00	See Figure 1 – step 1
Rinse	Running water	04:00	Rinse a let in a bath of running water See Figure 1 – step 1
Stain	Ponceau fuchsin	05:00	See Figure 1 – step 2
Rinse	Running water	01:00	See Figure 1 – step 2
Stain (mordant)	Phosphomolybdic acid	03:00	See Figure 1 – step 3
Stain	Light green	05:00	See Figure 1 – step 4
Rinse	1% acetic water	No	Quickly rinse the slide 2 successive baths See Figure 1 – step 5
Dehydrate	growing degree alcohols baths	No	until absolute alcohol
Dehydrate	Toluene or xylene	No	No
Mount	Toluene or Xylene based mounting medium	No	No

Figure 1. Kit Trichrome de Masson, light green variation staining steps

1 – 5: steps 1 to 5

## Expected results

**Nuclei:** black blue to brown

**Cytoplasm:** pink to red

**Erythrocytes:** bright pink

**Collagen and mucus:** green

**Elastic lamina:** pink

If observed results vary from those expected, please contact RAL Diagnostics technical service through your usual supplier for assistance.

## Performance

This medical device is state of the art. Its analytical performance, scientific validity and medical relevance are assessed in the CE marking review.

To ensure product performance, use clean and dry laboratory equipment.

The laboratory is responsible for notifying the manufacturer and state competent authority of any serious incident relating to the use of the medical device.

## User quality Control

Users are responsible for determining the appropriate quality control procedures for their laboratory and complying with applicable laboratory regulations.

RAL Diagnostics recommend quality control at reagents renewal and for the first staining cycle of each day. Slides stained for quality control purposes should be

checked to ensure that they are satisfactory for intended test (properly stained and free of precipitate).

Staining results for each cell type must also be compliant with this manual expected results

These quality control procedures should only be performed by qualified personnel.

## Other products

For more information contact your usual supplier.

## Recommendations, notes, and troubleshooting

### Products appearance

If the appearance of the products differs from the description above, do not use it and contact RAL Diagnostics technical service through your usual supplier for assistance.

### Procedures notes

To prevent products degradation, please comply with the storage and handling recommendations specified in this manual.

Staining times may vary according to the tissues section structure and the thickness of the section.

### Products stability

Every RAL Diagnostics product can be used until the expiry date indicated on, in its original packaging if it is still hermetically sealed.

### **Staining stability**

Staining quality and reproducibility depend on the correct use of the products. Staining conducted according to these recommendations will remain stable for several days. RAL Diagnostics recommends mounting the stained slides with a coverslip using a suitable mounting liquid and to store them in a light and dustproof container.

### **Instructions for cleaning and waste disposal**

All biological samples, effluents and used consumables should be considered potentially hazardous.



To avoid any risk, apply the following instructions: dispose of samples, effluents and consumables in accordance with laboratory standards and applicable national and local standards and regulations.

Chemical and biological waste must be collected and processed by specialized, registered companies.

## Table of symbols and abbreviations

Depending on the product, you may find the following symbols on the device or the packaging material.

GHS PICTOGRAMS	INTERPRETATION	SYMBOL	INTERPRETATION
	Explosive		Batch code
	Flammable		Serial number
	Oxidizer		Catalogue reference
	Compressed gas		Date of manufacture
	Corrosive		Use up to
	Toxic		Unique device identifier
	Harmful		Manufacturer
	Health Hazard		Importer
	Environmental Hazard		Entry distributing the medical advice in the region concerned
	No labelling applicable		CE marking device
			In vitro diagnostic medical device
			Authorised Representative in the European Community
			Authorised Representative in Switzerland
			Complies with UK guidelines
			Do not use if packaging is damaged
			Keep away from light
			Temperature limit: 15-25°C
			Temperature limit: 15-30°C
			Keep dry
			Box: handling upwards
			Fragile
			Sterilised by irradiation
			Single sterile barrier system with outer protective packaging
			Sterile and radiation-sterilised barrier suit
			Do not reuse
			Do not resterilize
			Contents sufficient for n tests
			Hazardous material contained
			Consult instructions for use
			Use
			After opening, use within XX months
			The product must not be used in conjunction with an automatic colouring machine
			Indicates a medical device that contains potentially carcinogenic, mutagenic or reprotoxic (CMR) substances, or substances classified as endocrine disruptors

## Bibliography

**CHEVREAU J., BELLOT J., CABANIER M-J.**, *Formulaire de Techniques Histologiques*, Maloine S.A. éd., 1977, p. 155-156.

**GANTER P., JOLLES G.**, *Histochimie normale et pathologique*, éd. Gauthier-Villars.,1970, vol. 2, p.1420.

## Change tracking

Date	Version	Changes
05/2022	IFU084A-RAL	IVDR (EU) 2017/746 compliance

