

Wide range of accessories

For optimal workflow

Wire basket with handle

The wire basket for convenient handling of autoclaving goods like solids, waste in bags, instruments, bottles, liquids or the like.

Accessory for	Material	Diameter	Height	Article number
Aluminum models	stainless steel	ø 23 cm	17 cm	8583010
Stainless steel models	stainless steel	ø 25 cm	22.5 cm	1300149



Instrument tray

The stable instrument tray with three feet and minimum water level indicator can also be used as a base plate for a wire basket or bottles. When purchasing a vertical autoclave, an instrument tray is already included in the scope of delivery.

Accessory for	Material	Diameter	Height	Article number
Aluminum models	aluminum	ø 22.5 cm	4 cm	1252203
Stainless steel models	stainless steel	ø 25 cm	4 cm	1300155



1252203

1300155

Tray Set with 3 trays

The tray set with 3 trays made of stainless steel is suitable for nutrient media, instruments, material testing, etc. When purchasing a horizontal autoclave, a tray set is already included in the scope of delivery.

Accessory for	Material		Width	Depth	Height	Article number
VacPro 12	stainless steel	Rack	16.8 cm	33.5 cm	16.6 cm	8500479
		Trays	13.8 cm	29.5 cm	1.5 cm	
		usable height	3 cm*	3 cm**	3 cm**	
VacPro 22	stainless steel	Rack	20 cm	40 cm	20 cm	8500481
		Trays	17 cm	35 cm	1.6 cm	
		usable height	3 cm*	4 cm**	4 cm**	
STE-45-T	stainless steel	Rack	25.9 cm	46 cm	26 cm	8500404
		Trays	22.6 cm	51 cm	1.7 cm	
		usable height	6.3 cm*	5.5 cm**	5.5 cm**	



*upper tray
**centre and bottom tray

Wire basket

The wire basket for convenient handling of autoclaving goods like solids, waste in bags, instruments, bottles, liquids or the like.

Accessory for	Material	Width	Depth	Height	Article number
VacPro 12	stainless steel	14.6 cm	30 cm	12.5 cm	8500401
VacPro 22	stainless steel	17.1 cm	38 cm	15.6 cm	8500403
STE-45-T	stainless steel	22.5 cm	51.1 cm	20.6 cm	8500406



Drip Tray

The drip tray is used to prevent chamber contamination by boiling over of liquids such as agar-agar, and as a base plate for bottles.

Accessory for	Material	Width	Depth	Height	Article number
VacPro 12	stainless steel	14.2 cm	29.8 cm	1.6 cm	1300765
VacPro 22	stainless steel	16.7 cm	37.8 cm	1.6 cm	1300767
STE-45-T	stainless steel	22.2 cm	51.1 cm	1.6 cm	1300768

