

Masques de protection 5 couches Norme KN95 équivalente FFP2



- en tissu non tissé
- couches filtrantes en coton et tissu non-tissé soufflé par fusion
- filtrent gouttelettes, poussières, bactéries, germes, fumée, et plus de 95% des nanoparticules > 0,3 µm
- pliables
- confortables, sûrs, respirants et hypoallergéniques
- clip de nez réglable, en métal
- 2 élastiques auriculaires
- taille unique
- en sachet de 5 masques

référence

Prix HT

TU3030 Masques 5 couches, les 5

Comparaison normes FFP2 et KN95

	NORME	rendement du filtre	Agent d'essai	Débit	Fuite vers l'intérieur	Résistance à l'air		dégagement de CO ₂
						inhalation	expiration	
FFP2	EN 149:2001+A1:2009 (norme européenne)	≥94%	NaCl et huile de paraffine	95 l/min	≤8%	≤500 Pa	≤300 Pa	≤1%
KN95	GB2626-2006 (norme chinoise)	≥95%	NaCl	85 l/min	≤8%	≤350 Pa	≤250 Pa	≤1%

Certificate of Compliance



No. 4Q200401G.DDE0D64

Certificate's Holder:

Dongguan Dacheng Electronic Technology Co., Ltd
No.11 Keyuan 1st Road, Tangxia Town, Dongguan City

Certification ECM Mark:



" uniquement pour les clients de Dasheng"

Product:
Model(s):

KN95 Mask
DC-KN95

Verification to:

Standard:
EN 149:2001+A1:2009
related to CE Directive(s):
R 2016/425 (Personal Protective Equipment)

Remark: This document has been issued on a voluntary basis and upon request of the manufacturer. It is our opinion that the technical documentation received from the manufacturer is satisfactory for the requirements of the ECM Certification Mark. The Conformity mark above can be affixed on the products accordingly to the ECM regulations about its release and its use.

Additional information and clarification about the Marking:



The manufacturer is responsible for the CE Making process, and if necessary, must refer to a Notified Body. This document has been issued on the basis of the regulation on ECM Voluntary Mark for the certification of products. RG01_ECM rev.3 available at: www.entecerma.it

Issuance date: 01 April 2020

Expiry date: 31 March 2025

Reviewer
Technical expert
Amanda Payne



Approver
ECM Service Director
Luca Bedonni

