



- jusqu'à +1000°C
- excellente résistance aux acides

Creusets filtrants



- avec filtre
- conformes DIN 12909

capacité	Ø sup x h	Ø base	référence	Prix HT
Filtre pores 6 µm				
8 ml	Ø 25 x h28 mm	Ø 20 mm	JP2760	
15 ml	Ø 30 x h35 mm	Ø 24 mm	JP2761	
25 ml	Ø 35 x h40 mm	Ø 28 mm	JP2762	
35 ml	Ø 40 x h43 mm	Ø 30 mm	JP2763	
50 ml	Ø 45 x h50 mm	Ø 36 mm	JP2764	
Filtre pores 7 µm				
8 ml	Ø 25 x h28 mm	Ø 20 mm	JP2770	
15 ml	Ø 30 x h35 mm	Ø 24 mm	JP2771	
25 ml	Ø 35 x h40 mm	Ø 28 mm	JP2772	
35 ml	Ø 40 x h43 mm	Ø 30 mm	JP2773	
50 ml	Ø 45 x h50 mm	Ø 36 mm	JP2774	
Filtre pores 8 µm				
8 ml	Ø 25 x h28 mm	Ø 20 mm	JP2780	
15 ml	Ø 30 x h35 mm	Ø 24 mm	JP2781	
25 ml	Ø 35 x h40 mm	Ø 28 mm	JP2782	
35 ml	Ø 40 x h43 mm	Ø 30 mm	JP2783	
50 ml	Ø 45 x h50 mm	Ø 36 mm	JP2784	

Entonnoirs filtrants de Büchner en porcelaine



- * conformes à la norme DIN 12905

capacité	Ø x h	longueur tige	pour papier filtre	référence	Prix HT
3 ml	Ø 18 x h47 mm	20 mm	Ø 19 mm	HG1270	
10 ml	Ø 27 x h60 mm	30 mm	Ø 29 mm	HG1271	
25 ml	Ø 40 x h76 mm	37 mm	Ø 42 mm	HG1272	
50 ml	* Ø 45 x h100 mm	41 mm	Ø 48 mm	HG1273	
75 ml	* Ø 55 x h129 mm	53 mm	Ø 50 mm	HG1274	
135 ml	* Ø 70 x h150 mm	75 mm	Ø 74 mm	HG1275	
290 ml	* Ø 90 x h165 mm	75 mm	Ø 95 mm	HG1276	
580 ml	* Ø 110 x h201 mm	90 mm	Ø 114 mm	HG1277	
795 ml	* Ø 125 x h232 mm	106 mm	Ø 130 mm	HG1278	
1250 ml	* Ø 150 x h235 mm	110 mm	Ø 155 mm	HG1279	
1900 ml	* Ø 185 x h260 mm	120 mm	Ø 190 mm	HG1280	
4300 ml	* Ø 240 x h330 mm	138 mm	Ø 250 mm	HG1281	
5800 ml	* Ø 270 x h300 mm	125 mm	Ø 279 mm	HG1282	
10600 ml	* Ø 320 x h350 mm	150 mm	Ø 330 mm	HG1283	

Entonnoirs filtrants de Büchner en porcelaine

- conformes à la norme DIN 12905



capacité	Ø x h	longueur tige	pour papier filtre	référence	Prix HT
10 ml	Ø 29 x h45 mm	28 mm	Ø 27 mm	HG1900	
25 ml	Ø 42 x h55 mm	35 mm	Ø 40 mm	HG1912	
35 ml	Ø 48 x h67 mm	43 mm	Ø 45 mm	HG1901	
70 ml	Ø 62 x h94 mm	64 mm	Ø 55 mm	HG1902	
120 ml	Ø 77 x h99 mm	64 mm	Ø 70 mm	HG1903	
240 ml	Ø 97 x h111 mm	71 mm	Ø 90 mm	HG1904	
400 ml	Ø 116 x h132 mm	83 mm	Ø 110 mm	HG1905	
600 ml	Ø 130 x h137 mm	85 mm	Ø 125 mm	HG1906	
1000 ml	Ø 156 x h155 mm	96 mm	Ø 150 mm	HG1907	
2000 ml	Ø 192 x h191 mm	117 mm	Ø 185 mm	HG1908	
4000 ml	Ø 248 x h207 mm	117 mm	Ø 240 mm	HG1909	
7000 ml	Ø 296 x h235 mm	130 mm	Ø 270 mm	HG1910	
10000 ml	Ø 334 x h253 mm	133 mm	Ø 320 mm	HG1911	

Entonnoirs filtrants en porcelaine



- conformes DIN 12905

capacité	Ø haut. / base x h	Ø tige	pour papier-filtre	référence	Prix HT
50 ml	30 / 8 x 45 mm	7 mm	Ø 9 mm	HG1951	
100 ml	50 / 16 x 74 mm	11 mm	Ø 15 mm	HG1952	
150 ml	70 / 25 x 102 mm	16 mm	Ø 27 mm	HG1953	