


Tuyaux caoutchouc

- dureté Shore : 50 A
- tuyaux rouges



Ø int. mm	Ø ext. mm	paroi mm	réf.	Prix HT le mètre	réf.	Prix HT les 25 mètres
Ø 2	Ø 4	1,0	ND1001	-	-	-
Ø 3	Ø 5	1,0	ND1002	-	-	-
Ø 3	Ø 6	1,5	ND1003	-	-	-
Ø 4	Ø 6	1,0	ND1004	-	-	-
Ø 4	Ø 7	1,5	ND1005	-	-	-
Ø 4	Ø 8	2,0	ND1006	-	ND100625	-
Ø 5	Ø 8	1,5	ND1007	-	ND100725	-
Ø 5	Ø 9	2,0	ND1008	-	-	-
Ø 6	Ø 9	1,5	ND1009	-	ND100925	-
Ø 6	Ø 10	2,0	ND1010	-	ND101025	-
Ø 7	Ø 10	1,5	ND1011	-	-	-
Ø 7	Ø 11	2,0	ND1012	-	ND101225	-
Ø 8	Ø 12	2,0	ND1013	-	-	-
Ø 9	Ø 13	2,0	ND1014	-	ND101425	-
Ø 10	Ø 14	2,0	ND1015	-	ND101525	-
Ø 12	Ø 17	2,5	ND1016	-	ND101625	-
Ø 13	Ø 19	3,0	ND1017	-	-	-
Ø 14	Ø 20	3,0	ND1018	-	-	-
Ø 15	Ø 21	3,0	ND1019	-	-	-
Ø 16	Ø 22	3,0	ND1020	-	-	-
Ø 18	Ø 24	3,0	ND1021	-	-	-
Ø 20	Ø 28	4,0	ND1022	-	-	-
Ø 25	Ø 33	4,0	ND1023	-	-	-


Tuyaux à vide en caoutchouc

Ø int. mm	Ø ext. mm	paroi mm	réf.	Prix HT le mètre	réf.	Prix HT les 25 mètres
Ø 3	Ø 9	3	ND1051	-	-	-
Ø 4	Ø 12	4	ND1052	-	ND105225	-
Ø 5	Ø 13	4	ND1053	-	-	-
Ø 5	Ø 15	5	ND1054	-	-	-
Ø 6	Ø 12	3	ND1055	-	-	-
Ø 6	Ø 14	4	ND1056	-	-	-
Ø 6	Ø 16	5	ND1057	-	-	-
Ø 6	Ø 18	6	ND1058	-	-	-
Ø 7	Ø 15	4	ND1059	-	-	-
Ø 7	Ø 17	5	ND1060	-	-	-
Ø 8	Ø 16	4	ND1061	-	-	-
Ø 8	Ø 18	5	ND1062	-	ND106225	-
Ø 8	Ø 20	6	ND1063	-	-	-
Ø 8	Ø 24	8	ND1064	-	-	-
Ø 10	Ø 16	3	ND1065	-	-	-
Ø 10	Ø 20	5	ND1066	-	-	-
Ø 10	Ø 26	8	ND1067	-	-	-
Ø 10	Ø 30	10	ND1068	-	-	-
Ø 12	Ø 22	5	ND1069	-	-	-
Ø 12	Ø 28	8	ND1070	-	-	-
Ø 15	Ø 31	8	ND1071	-	-	-

Tuyaux de laboratoire

Tuyaux en PVC

- tuyaux en PVC transparents
- température : -20°C à +60°C
- le mètre



Ø int. mm	Ø ext. mm	paroi mm	réf.	Prix HT le mètre
Ø 1	Ø 1,8	0,40	ND1102	-
Ø 1,5	Ø 2,1	0,30	ND1103	-
Ø 2	Ø 3,0	0,50	ND1105	-
Ø 2	Ø 4,0	1,00	ND1106	-
Ø 3	Ø 4,0	0,50	ND1107	-
Ø 3	Ø 5,0	1,00	ND1108	-
Ø 3	Ø 6,0	1,50	ND1109	-
Ø 4	Ø 5,0	0,50	ND1110	-
Ø 4	Ø 6,0	1,00	ND1111	-
Ø 4	Ø 7,0	1,50	ND1112	-
Ø 4	Ø 8,0	2,00	ND1113	-
Ø 5	Ø 7,0	1,00	ND1114	-
Ø 5	Ø 8,0	1,50	ND1115	-
Ø 5	Ø 9,0	2,00	ND1116	-
Ø 6	Ø 8,0	1,00	ND1117	-
Ø 6	Ø 9,0	1,50	ND1118	-
Ø 8	Ø 10,0	1,00	ND1122	-
Ø 6	Ø 10,0	2,0	ND1119	-
Ø 7	Ø 10,0	1,5	ND1120	-
Ø 7	Ø 11,0	2,0	ND1121	-
Ø 8	Ø 11,0	1,5	ND1123	-
Ø 8	Ø 12,0	2,0	ND1124	-
Ø 9	Ø 12,0	1,5	ND1125	-
Ø 9	Ø 13,0	2,0	ND1126	-
Ø 10	Ø 13,0	1,5	ND1127	-
Ø 10	Ø 14,0	2,0	ND1128	-
Ø 10	Ø 16,0	3,0	ND1129	-
Ø 12	Ø 15,0	1,5	ND1130	-
Ø 12	Ø 16,0	2,0	ND1131	-
Ø 12	Ø 18,0	3,0	ND1132	-
Ø 13	Ø 17,0	2,0	ND1133	-
Ø 13	Ø 19,0	3,0	ND1134	-
Ø 15	Ø 19,0	2,0	ND1135	-
Ø 16	Ø 20,0	2,0	ND1136	-
Ø 18	Ø 23,0	2,5	ND1146	-
Ø 19	Ø 27,0	4,0	ND1139	-
Ø 20	Ø 26,0	3,0	ND1140	-
Ø 22	Ø 28,0	3,0	ND1141	-
Ø 25	Ø 34,0	4,0	ND1142	-
Ø 30	Ø 39,0	4,5	ND1143	-
Ø 35	Ø 45,0	5,0	ND1144	-
Ø 40	Ø 50,0	5,0	ND1145	-

Ø int. mm	Ø ext. mm	paroi mm	réf.	Prix HT les 50 mètres
Ø 8	Ø 10,0	1,00	ND112250	-
Ø 8	Ø 12,0	2,0	ND112450	-
Ø 9	Ø 12,0	1,5	ND112550	-
Ø 10	Ø 14,0	2,0	ND112850	-
Ø 10	Ø 16,0	3,0	ND112950	-
Ø 15	Ø 19,0	2,0	ND113550	-

Ø int. mm	Ø ext. mm	paroi mm	réf.	Prix HT les 100 mètres
Ø 2	Ø 4,0	1,00	ND1106100	-
Ø 3	Ø 4,0	0,50	ND1107100	-
Ø 3	Ø 5,0	1,00	ND1108100	-
Ø 3	Ø 6,0	1,50	ND1109100	-
Ø 4	Ø 7,0	1,50	ND1112100	-
Ø 5	Ø 8,0	1,50	ND1115100	-
Ø 6	Ø 8,0	1,00	ND1117100	-